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
User Control Variables			
		Minimum	Maximum
CV 4	Image acq. delay (sec)	0	100
CV 16	Uni-polar echoes (0=off, 1=on)	0	1

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New Protocol
E-SWAN

Processing Thresholds

- Set sliders to select range of MR Units to be processed.

Pixels within the selected range of values are shown in green.

Filter Width

Width of the low-pass frequency used to separate phase shift caused by background.

The value must be one positive integer.

Skip Echo Indices

Ignore the echo when do calculation. Echo index starts from 1.

The syntax should be something like "1 3 9 ---", which means that Echo No.1, No.3 and No.9 will be skipped when do calculation.

Confidence Level

The confidence level parameter is small in areas corresponding to a high "confidence" and large in areas corresponding to a low "confidence".

Confidence Level only applies if more than two images at one location are provided.

Processing Steps:

- Step 1/5: Processing Thresholds (Noise slider from 1 to 32759)
- Step 2/5: Filter Width (Radio buttons: 4, 8, 16, 32, 64, Custom)
- Step 3/5: Skip Echo Indices (Text input: ---)
- Step 4/5: Confidence Level (Radio buttons: 0.010, 0.05, 0.10, 0.30, 0.80, Custom)

ESWAN Functool后处理 (Magnitude, Phase) :

Functool 5.4.07

New Protocol
E-SWAN

Final Settings

Threshold	Noise 1
	Upper 32759
filter-width	32
skip-channel	---
conf-level	0.05

Compute

Close Back 5/5 Next

Time Graph


Phase

Color Ramps

- Display Normal
- Save View
- SaveView + AddToReport
- Create Annotation
- Clear Region
- Split View

Rainbow
Hot Iron
Puh thalam
3 colors
Gray levels
Inv. Rainbow

ESWAN Functool后处理 (R2STAR, T2STAR):

Confidence Level

The confidence level parameter is small in areas corresponding to a high "confidence" and large in areas corresponding to a low "confidence".

Confidence Level only applies if more than two images at one location are provided.

Confidence level: 0.010 0.05 0.10 0.30 0.80 Custom

Final Settings

Threshold	Noise 1
Upper 32759	
filter-width	32
skip-channel	---
conf-level	0.05

Buttons: Compute, Back, 5/5, Next

ESWAN Functool 批量保存图像:

The screenshot displays the GE Functool 5.4.07 software interface. The main window shows an axial MRI brain scan with two regions of interest (ROIs) highlighted in purple and green. The interface includes a toolbar on the left with various navigation and analysis tools, and a data panel on the right showing ROI statistics.

ROI Statistics (Left Panel):

ROI	min	Avg.	Dev.
1	168	3186.8	352.04
2	132	3131.5	127.75

ROI Statistics (Right Panel):

ROI	min	Avg.	Dev.
1	158	3156.8	352.04
2	132	3131.5	127.75

File / Save Functional Map(s) Dialog:

Select functional map(s):

- Visible maps
 - Left
 - Right
 - All
- Select from available maps...

Select single or multiple location(s):

- Single location
- Multiple locations

File / Save Format Dialog:

Select format:

- Send to Film Composer (F1)
- Save as SCREENSAVE image (S)
- SCREENSAVE + Add To Report (D)
- Save as processed images (F) (PROCESS / REFORMAT / SCPT)

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